




Erisma[®] - LP MIS

Pedicle Screw System




3-ZONE THREAD

Initiation Zone with aggressive tip and dual-lead thread for ease of insertion. Cancellous and Cortical Zones for ideal bone purchase.



ROBUST TABS

Tabs remain reliably intact until easily removed with the Tab Breaker.



RELIABLE & PRECISE ROD PLACEMENT

Low-profile inserter securely retains rod for controlled placement in-situ.



PATENTED PERSUADER DESIGN

Controlled rod persuasion with integrated set screw insertion.

SIZES*

Erisma® - LP MIS



POLYAXIAL SCREWS

4.5mm x 20 - 45mm	8.5mm x 25 - 100mm
5.5mm x 20 - 55mm	9.5mm x 25 - 100mm
6.5mm x 20 - 60mm	10.5mm x 25 - 100mm
7.5mm x 25 - 100mm	



SET SCREW



PRE-BENT RODS



STRAIGHT RODS

Prior to using any Clariance device, please review the instructions for use and surgical technique for a complete listing of indications, contraindications, warnings, precautions, potential adverse events, and directions for use.

Clariance, 18 Rue Robespierre 62217 Beaurains France • T : +33 3 2116 12 15 F : +33 3 2115 50 73 • contact@clariance-spine.com
 Clariance, Inc 4809 N. Ravenswood Ave Suite 119 Chicago, IL 60640 • T : 773-868-7041 F : 773-868-7043 • contact@clariance-spine.us



*Product availability may vary upon markets.
 09/2019 Clariance, Inc All rights reserved. Clariance is a registered trademark licensed by Clariance, Inc. Idys is a trademark licensed by Clariance Inc - 29FLY01